




The *Stable* Isotope Company

EA Autosamplers





Standard 66 place
Semi-open carousel design, samples
individually purged to exclude
atmosphere.




Zero Blank 66 place
Sealed carousel design to exclude
atmosphere from samples
completely.

The key features of the SerCon EA Autosamplers :

- 66 place autosampler for unattended operation.
- 130 place multilayer option available.
- High quality all metal construction.
- Pneumatic operation (3 bar).
- Replaceable piston and sleeve assembly on Standard autosampler.
- Internal shutoff valve on Zero Blank autosampler protects combustion reactor.
- Optional support plate.
- Compatible with SerCon, Europa, CarloErba and other EAs.
- Zero Blank autosampler is suitable for use on TC-EA pyrolysis systems.



Standard Autosampler



Zero Blank Autosampler

