




The *Stable* Isotope Company

## High Temp EA-Liquid



SerCon present the High Temp EA, the most robust 1450°C thermal decomposition elemental analyzer now with liquid injection capability


Providing High Precision  $^{18}\text{O}$  and  $^2\text{H}$  sample analysis robustly

Designed for utmost reliability from 600°C to 1450°C, analysis without failure

Designed to operate with the revolutionary SerCon 20-20 IRMS and SerCon Elemental Analysers

Compatible with all elemental analyzer and IRMS combinations

The High Temp EA provides the ultimate performance



Actual Screen Graphic from SerCon Callisto Software

